

# RRDE-3A Rotator

## For Electrochemical Measurements

### Feature of RRDE-3A Rotator

- Can be used for both **Rotating Disk Electrode** and **Rotating Ring Disk Electrodes** measurements
- Remote and manual control modes
- Motor rotating speed range: 90-8,000 rpm
- Equipped with spin coating function
- Equipped with gas purge function
- Easy for electrodes and cell vial exchange
- Low mechanical noise
- High accuracy of rotating speed
- Custom-made electrodes available
- Compact design & Easy operation



## RRDE-3A Rotating Ring Disk Electrode Rotator



Spin- Coating State



RDE & RRDE electrodes for RRDE-3A Rotator